









INJECTOR PORT LINERS – AGILENT GCs

All injection port liners are polysiloxane deactivated and ready for use. They are manufactured to exacting tolerances and come in various configurations for many popular GCs. Guaranteed inert for phenols, organic acids and bases.

CONFIGURATION	DIMENSIONS	OEM PART #	QTY	QUADREX CAT. #	PRICE
 Splitless	2mm I.D. X 79mm	5181-8818 18740-80220	1	DGL-208111	\$25.36
			10	DGL-208111-10	\$226.00
 Split/splitless	4mm I.D. X 79mm	210-3003	1	DGL-208112	\$15.44
			10	DGL-208112-10	\$109.15
 Split/splitless Glass wool packed	4mm I.D. X 79mm	19251-60540	1	DGL-208113	\$15.44
			10	DGL-208113-10	\$138.00
 Cup Splitter	4mm I.D. X 79mm	18740-80190	1	DGL-208114	\$48.51
			10	DGL-208114-10	\$411.23
 Split/splitless Cup, packed with 10% OV-1 on Chromosorb WHP	4mm I.D. X 79mm	187-60840	1	DGL-208115	\$54.02
			10	DGL-208115-10	\$475.18
 Dual tapered	4mm I.D. X 79mm	5181-3315	1	DGL-208116	\$28.67
			10	DGL-208116-10	\$243.65
 Single taper	4mm I.D. X 79mm	5181-3316	1	DGL-208119	\$19.85
			10	DGL-208119-10	\$179.71
 Single taper with glass wool	4mm I.D. X 79mm	5062-3587	1	DGL-208120	\$28.67
			10	DGL-208120-10	\$241.83