






## INJECTOR PORT LINERS – PERKIN-ELMER GCs

All injection port liners are polysiloxane deactivated and ready for use (other deactivation available, inquiry for details). They are manufactured to exacting tolerances and come in various configurations for many popular GCs. Guaranteed inert for phenols, organic acids and bases.

CONFIGURATION	DIMENSIONS	OEM PART #	QTY	QUADREX CAT. #	PRICE
 'Dimple' Split	100mm	0330-5181	1	<b>DGL-208212</b>	<b>\$26.46</b>
			10	<b>DGL-208212-10</b>	<b>\$220.50</b>
 Splitless	100mm	0330-5180	1	<b>DGL-208214</b>	<b>\$41.90</b>
			10	<b>DGL-208214-10</b>	<b>\$359.42</b>
 Split/Autosystem	4mm I.D. X 92mm	N610-1052	1	<b>DGL-208232</b>	<b>\$24.26</b>
			10	<b>DGL-208232-10</b>	<b>\$197.35</b>
 Split/Autosystem w/ Deactivated Glass Wool	4mm I.D. X 92mm		1	<b>DGL-208233</b>	<b>\$27.56</b>
			10	<b>DGL-208233-10</b>	<b>\$246.96</b>
 Split/Autosystem	2mm I.D. X 92mm	N612-1372	1	<b>DGL-208234</b>	<b>\$24.26</b>
			10	<b>DGL-208234-10</b>	<b>\$185.22</b>