



INJECTOR PORT LINERS – THERMO-FINNIGAN GCs

All injection port liners are polysiloxane deactivated and ready for use (other deactivation available, inquiry for details). They are manufactured to exacting tolerances and come in various configurations for many popular GCs. Guaranteed inert for phenols, organic acids and bases.

CONFIGURATION		DIMENSIONS	OEM PART #	QTY	QUADREX CAT. #	PRICE
	Split	3mm I.D. x 105mm	453 20031	1	DGL-208611	\$24.26
				10	DGL-208611-10	\$209.48
	Split	5mm I.D. x 105mm	453 20030	1	DGL-208612	\$22.05
				10	DGL-208612-10	\$185.22
	Splitless	3mm I.D. x 105mm	453 20032	1	DGL-208615	\$36.38
				10	DGL-208615-10	\$319.73
	Splitless	5mm I.D. x 105mm	453 20033	1	DGL-208616	\$35.28
				10	DGL-208616-10	\$303.19