







## INJECTOR PORT LINERS – VARIAN GCs

All injection port liners are polysiloxane deactivated and ready for use (other deactivation available, inquiry for details). They are manufactured to exacting tolerances and come in various configurations for many popular GCs. Guaranteed inert for phenols, organic acids and bases.

CONFIGURATION		DIMENSIONS	OEM PART #	QTY	QUADREX CAT. #	PRICE
	Open Split	4mm I.D. X 72mm	16-000830-00	1	<b>DGL-208311</b>	<b>\$24.26</b>
				10	<b>DGL-208311-10</b>	<b>\$190.73</b>
	Frit Split	4mm I.D. X 72mm	01-900109-03 16-000830-01	1	<b>DGL-208312</b>	<b>\$39.69</b>
				10	<b>DGL-208312-10</b>	<b>\$329.65</b>
	Open Split w/ Glass Wool packed	4mm I.D. X 72mm	01-900109-01	1	<b>DGL-208313</b>	<b>\$28.67</b>
				10	<b>DGL-208313-10</b>	<b>\$231.53</b>
	Splitless Borosilicate	2mm I.D. X 74mm	01-900109-05 03-949437-00	1	<b>DGL-208314</b>	<b>\$27.56</b>
				10	<b>DGL-208314-10</b>	<b>\$213.89</b>
	SPI	0.5mm I.D. X 54mm	01-900109-06 03-918332-01	1	<b>DGL-208316</b>	<b>\$40.79</b>
				10	<b>DGL-208316-10</b>	<b>\$300.98</b>
	SPI	0.8mm I.D. X 54mm	01-900109-07 03-918332-02	1	<b>DGL-208318</b>	<b>\$40.79</b>
				10	<b>DGL-208318-10</b>	<b>\$300.98</b>