




GC CAPILLARY COLUMN

APPLICATIONS: **Ultra Alloy™** **Stainless Steel**

Sesame Oil

COLUMN: UAC-TRG, Triglyceride Phase
30M. x 0.25mm I.D. x 0.15µm film
Cat. No.: UAC-TRG-30-0.15F

Temperature: 340° (1°/min.) - 360°C



QUADREX



PO Box 3881, Woodbridge, CT 06525 • (203) 393-3112 • FAX (203) 393-0391
TOLLFREE: 800-275-7033 • E-MAIL: tech@quadrexcorp.com • WEB: www.quadrexcorp.com